


95 % **BSI sCMOS** High Quantum Efficiency

Scientific CMOS Camera

Dhyana 95






Backside-illuminated sCMOS technology Opening a new era of high sensitivity imaging applications!

The Dhyana 95 is a highly sensitive scientific camera developed around back-illuminated sCMOS technology, not only with ultra-high quantum efficiency comparable to that of an EMCCD, but faster frame rate and excellent signal-to-noise performance, as well as a 2" field of view, 200-1100nm spectral response and dynamic range of 11um pixels.

These innovative features open a new era to disciplines such as: life sciences; spectral analysis; astronomical observation and other cutting-edge research areas.







11µm x 11µm Pixel BSI sCMOS sensor 1.6e-(Median) Low Read Noise



90,000e-Full Well Capacity




-15°C Cooling Low Dark Current



USB3.0 High Frame Speed

Ultra-high Quantum Efficiency


Quantum efficiency refers to the ratio of the average number of photoelectrons produced per unit time at a particular wavelength to the number of incident photons. It is an important parameter to describe the photoelectric conversion capability of optoelectronic devices and is one of the important indicators to measure camera sensitivity.




The Dhyana95 uses backside-illuminated sCMOS thinned chip technology to avoid light interference from the wiring layer, thereby increasing the pixel area and improving the photoelectric conversion rate, as shown in the figure: at 560 nm, the quantum efficiency is 95%, which a front illuminated sCMOS camera cannot match, and is comparable to the best EMCCD camera.

Excellent Signal to Noise Ratio

Signal-to-noise ratio is the ratio of detectable signal to background noise. The higher the signal-to-noise ratio, the more obvious the recognizable signal is.



The Dhyana95 read noise is only 1.6 electrons(Median), therefore the signal to noise ratio is significantly better than other sCMOS cameras, and when the incident photons are>3, there is a better performance than EMCCD based cameras.




2" Large Field of View

The 2" array can not only adapt to more optical interfaces and deliver a greater field of view, but also results in fewer lens switches to find the area of interest on the sample.

- 2" Tucsen Dhyana 95
- 1.2" sCMOS Camera
- 2/3" EMCCD Camera

Challenge the Most Extreme Applications

The Dhyana95's significant breakthrough in quantum efficiency means that sCMOS cameras have the potential to challenge much more extreme applications, with lower excitation energy, lower dye concentrations, and weaker photon signals, and still provide rich image information.



Tucsen Dhyana 95

The 97% QE EMCCD camera

The 82% QE sCMOS camera


Under very weak laser illumination (< 0.002 kW cm⁻²), Dhyana 95 camera showed comparable fluorescent intensity to the 97%QE EMCCD camera with 100 EM gain. On the other hand, Dhyana 95 had similar background fluorescence noise when comparing to the 82%QE sCMOS camera, which has 1.0e-low read noise.

Customer Applications

"Extensive tests and comparisons with other top-of-the-line EMCCD and sCMOS cameras have been performed in our laboratory. We found Dhyana 95, the new sCMOS camera stood its own ground remarkably well and offered satisfactory performance across the board". —Professor Ning Fang at Georgia State University

STORM Imaging

"The beta tubulin of Hela cell was immunostained with mouse anti-tubulin primary antibody and Alexa fluor 647 conjugated Donkey anti mouse secondary antibody. The ROI is about 50μ m x 50μ m. For Epifluorescence images the exposure time is 100 ms giving imaging rate of 10 fps. It is 50 fps for STORM imaging and was taken over 40000 frames to reconstruct the STORM images."



The localization precision from Dhyana 95 was calculated to be 6.9 nm, an increase of 8% and 35% in localization precision compared to the 97%QE EMCCD camera (7.5 nm) and the 82%QE sCMOS camera (9.3 nm), respectively.

a. Tucsen Dhyana 95 b. the 97% QE EMCCD camera c. the 82% QE sCMOS camera

The FRC resolution of the Tucsen Dhyana 95 camera (40nm) was comparable to that (39nm) of the 97%QE EMCCD camera and slightly better than that (44nm) of the 82%QE sCMOS camera.

a. Tucsen Dhyana 95

- b. the 97% QE EMCCD camera
- c. the 82% QE sCMOS camera

Technical Features

Model	Dhyana 95
Sensor size	2"
Sensor model	Backside-illuminated sCMOS
Color/monochrome	Monochrome
Quantum efficency	79%@240nm, 95%@560nm, 60%@800nm
Effective no.of pixels	2048(H) x 2048(V)
Pixel size	11µm x 11µm
Effective area	22.528mm x 22.528mm
Full well capacity	90,000e-
Frame rate (full resolution)	24fps@16bit
Read noise	1.6e-(Median)
Shutter type	Rolling Shutter
Exposure time	0.021ms-10s
Cooling method	Peltier cooling
Cooling temperature	Forced air (Ambient at +20°C): -15°C
Dark current	3 electrons / pixel / s (0°C) (typ.) 1.5 electrons / pixel / s (-15°C) (typ.)
DSNU	0.2e-
Sub-array	Available
Sub-array External trigger mode	Available Standard / Synchronous / Global trigger
External trigger mode	Standard / Synchronous / Global trigger
External trigger mode Trigger delay function	Standard / Synchronous / Global trigger 0-10,000s
External trigger mode Trigger delay function Signal output ports	Standard / Synchronous / Global trigger 0-10,000s 3 (Exposure / Global / Readout signal)
External trigger mode Trigger delay function Signal output ports External trigger connector	Standard / Synchronous / Global trigger 0-10,000s 3 (Exposure / Global / Readout signal) SMA
External trigger mode Trigger delay function Signal output ports External trigger connector Digital interface	Standard / Synchronous / Global trigger 0-10,000s 3 (Exposure / Global / Readout signal) SMA USB3.0
External trigger mode Trigger delay function Signal output ports External trigger connector Digital interface SDK	Standard / Synchronous / Global trigger 0-10,000s 3 (Exposure / Global / Readout signal) SMA USB3.0 Support
External trigger mode Trigger delay function Signal output ports External trigger connector Digital interface SDK Bit depth	Standard / Synchronous / Global trigger 0-10,000s 3 (Exposure / Global / Readout signal) SMA USB3.0 Support 16 bit
External trigger mode Trigger delay function Signal output ports External trigger connector Digital interface SDK Bit depth Lens mount	Standard / Synchronous / Global trigger 0-10,000s 3 (Exposure / Global / Readout signal) SMA USB3.0 Support 16 bit T or C-mount
External trigger mode Trigger delay function Signal output ports External trigger connector Digital interface SDK Bit depth Lens mount Power supply	Standard / Synchronous / Global trigger 0-10,000s 3 (Exposure / Global / Readout signal) SMA USB3.0 Support 16 bit T or C-mount 12V / 8A
External trigger mode Trigger delay function Signal output ports External trigger connector Digital interface SDK Bit depth Lens mount Power supply Power consumption	Standard / Synchronous / Global trigger 0-10,000s 3 (Exposure / Global / Readout signal) SMA USB3.0 Support 16 bit T or C-mount 12V / 8A 50W
External trigger mode Trigger delay function Signal output ports External trigger connector Digital interface Bit depth Lens mount Power supply Power consumption Camera size	Standard / Synchronous / Global trigger0-10,000s3 (Exposure / Global / Readout signal)SMAUSB3.0Support16 bitT or C-mount12V / 8A50W120mm x 119mm x 121mm
External trigger mode Trigger delay function Signal output ports External trigger connector Digital interface SDK Bit depth Lens mount Power supply Power consumption Camera size PC software	Standard / Synchronous / Global trigger 0-10,000s 3 (Exposure / Global / Readout signal) SMA USB3.0 Support 16 bit T or C-mount 12V / 8A 50W 120mm x 119mm x 121mm Mosaic/LabVIEW/Matlab/Micromanager Windows/Linux is supported ,

Application

- Super-resolution microscopy
- Real-time confocal microscopy
- Gene sequencing
- Live-cell imaging
- Single molecule detection
- Astronomical observation
- TIRF
- FRET

Dimensions


118.74 70.71

M42-P0.75

4*M4DP3


0

FØ100.00



0

်ဝ



Tucsen Photonics Co., Ltd.

.00

Address: 6F NO.1 building Caimao Zone, 756# Qi'an Road, Gaishan Town, Cangshan Area, Fuzhou, Fujian, P. R. CHINA Tel: + 86-591-88194580 Website: www.tucsen.com Email: support@tucsen.com Europe/North America 0EM/0DM Support:

Address: 1273 Vine St, Denver, CO 80206, United States Tel: +1-520-203-2643 E-mail: mikeblake@tucsen.com